

## Wide range of application



**Scene 1**  
Infusion status control



**Scene 3**  
Automatic calculation of IV drip




**Scene 2**  
Infusion status monitoring



**Scene 4**  
Adjustment support



Responding to on-site needs,  
providing safety and peace of mind for users



### Specifications

Application	Infusion set *Details to be determined	External dimensions	148(W)×149(H)×41(D)mm
Product name	Infusion multi-monitoring system(Tentative)	Weight	310g
Measurement range	Infusion volume 0 - 9999mL Infusion time 0 - 9999分 Measurement accuracy *unverified Measurement range	Communication module	802.11 b/g/n 2.4-2.5GHz(technical standards compliance)
Battery	Lithium-ion battery, About 5 hours	Application	Android application rowser-based client application Browser-Based Server application Wi-Fi settings(Captive Portal) Initial setup PC software
Notification in case of abnormalities	No start operation performed / No operation Performed after completion / Drip rate threshold exceeded / Drip not detected Overtime / Low battery / Data error/ Loss of Communication		



Contact **P&M Co.,Ltd.**

11. Monden-machi Industrial Park, Aizuwakamatu-city,  
FUKUSHIMA, 965-0845, JAPAN

Email: sales@pandm-aizu.com

